# wiroMEDI LOOP

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# wiroMEDI LOOP

## MEDICAL FACE MASKS TYPE IIR

Polypropylene spun-bonded fabric | 3-ply | Ear loops | Nosepiece

PREMIUM

# Product data sheet

## wiroMEDI LOOP

### Manufacturer


Company + Address	WIROS Wilfried Rosbach GmbH, Karl-Arnold-Str. 5, 47877 Willich, Germany
SRN	DE-MF-000004956
Trademark	care & serve®

### Authorisation

CH REP	Swiss AR Services AG, Industriestrasse 47, 6300 Zug, Switzerland
--------	--

### Product description

Category	PREMIUM
Product description	care & serve® PREMIUM   wiroMEDI LOOP   Medical face masks type IIR Polypropylene spun-bonded fabric   3-ply   Ear loops   Nosepiece
Product features	Bacterial filtration efficiency (BFE) > 99% 3-ply ear loops nosepiece

### Material

Material	Non-woven polypropylene (3-ply)
Material composition	internal
Material weight	Outer layer: $\approx 28 \text{ g/m}^2 (\pm 2 \text{ g/m}^2)$ Filter medium: $\approx 25 \text{ g/m}^2 (\pm 2 \text{ g/m}^2)$ Inner layer: $\approx 23 \text{ g/m}^2 (\pm 2 \text{ g/m}^2)$

### Application

Intended use	Disposable medical face mask for covering the mouth and nose to protect patients against infectious germs, minimise the direct transfer of infectious germs between personnel and patients and protect the wearer against contaminated splashes of liquid. Medical face masks can also be worn by patients and other individuals to reduce the risk of spreading infection, particularly in epidemic or pandemic situations.
Contact with food produce	Short-term indirect contact with foodstuffs is permitted.
Scope of application	Face masks for protection against the transfer of infectious germs during medical applications and procedures as well as in epidemic or pandemic situations.
Application group	Use in the medical field is restricted to appropriately trained professionals. Also suitable for laypersons in an epidemic and/or pandemic situation.
Patient group	Face masks are suitable for all patient groups provided the mask fits well while worn and is not too large or small.
Indication	Face masks for protection against the transfer of infectious germs during medical applications and procedures as well as in epidemic or pandemic situations.
Contraindication	No known contraindications.
Notes on use	Disposable product Not sterile
Limitation	The nose clip of the masks (face masks) is magnetic and not suitable in combination with highly magnetic devices, e.g. MRT devices. Not suitable for packing or storing food.
Note	All serious incidents directly related to this product must be reported to the manufacturer and the competent authorities.

# Product data sheet

## wiroMEDI LOOP

Application areas	Hospital   Old people's home   inc. specialist circles Medical practices Industry   Pharmaceutical industry Laboratory Food processing Cleaning   Hygiene Cosmetics   Wellness   Beauty industry	      
Supplementary material info	Latex-free	
Health advice	<p>If skin irritation or allergic reactions occur, discontinue use immediately and seek medical advice.</p> <p>As far as is known according to the current state of the art, the products do not contain any toxic, carcinogenic or mutagenic substances when used in accordance with their intended use.</p> <p>These products do not contain latex components and are considered hypoallergenic. However, they can cause irritation in sensitive skin. Long-term wear may restrict skin ventilation, which can cause skin irritation.</p>	

### Standards

CE marking	Class I medical device in accordance with EU 2017/745	 								
SRN	DE-MF-000004956									
Basis UDI DI	4051642-wiroMEDILoop-JD									
	<p>CE marking</p> <p>The CE mark guarantees the free movement of products and goods within the European Economic Community. All products with the CE mark comply with the essential requirements of the European regulations specified.</p>									
EN 14683:2019+AC:2019	<p>Medical face masks – requirements and test methods</p> <table border="1"> <thead> <tr> <th>Test</th> <th>Result</th> </tr> </thead> <tbody> <tr> <td>Bacterial filtration efficiency (BFE), (%)</td> <td>&gt; 99%</td> </tr> <tr> <td>Pressure difference (Pa/cm<sup>2</sup>)</td> <td>&lt; 60</td> </tr> <tr> <td>Microbiological purity (KBE/g)</td> <td>≤ 30</td> </tr> </tbody> </table> <p>Information</p> <p>User information</p> <p>When breathing, coughing and sneezing, smaller or larger quantities of respiratory droplet are released from the mucous membrane in the mouth and nose. Most droplets range from 0.5 µm to 12 µm in diameter. Larger droplets in particular can contain microorganisms from aerosolised particles. The droplets can then spread by air to sensitive places, such as an open surgical wound or sterile instruments.</p> <p>Medical masks, which are intended for use in operating theatres and healthcare facilities with similar requirements, are designed to protect the entire working environment. This standard describes two types of medical masks and their respective levels of protection. A Type I medical mask is used around patients to reduce the risk of spreading infection, particularly in epidemic or pandemic situations. A Type II mask is primarily intended for use by medical staff in an operating theatre or other medical facility with similar requirements.</p> <p>A special situation, which is also handled in the European Medical Device Regulation, is one in which the wearer requires spray protection against potentially contaminated liquids. If the mask is intended to protect the wearer against infectious germs (bacteria, viruses and fungi) when used correctly, respiratory protective equipment must be attached in accordance with Directive (89/686/EEC) and Regulation (EU) 2016/425 on personal protective equipment (PPE). The performance requirements of respiratory protective equipment are covered in EN 149.</p> <p>The effectiveness of a mask is dependent on a number of factors, such as filtration efficiency and how the mask sits on the wearer's face. Different designs suit different areas of application and so it is important to carefully choose a mask that produces the desired result.</p> <p>The filtration efficiency of mask materials can differ depending on the filter media. How the mask fits varies considerably between those who secure it behind their ears using ear loops and those who secure it using straps over the head and a nose piece that can be adjusted to fit the shape of the face. Mask performance is typically determined by in-vitro tests of the materials from which the mask is made. However, it is important to carefully consider how the mask fits when choosing a mask for a specific purpose. Another factor to consider is the mask's ability to absorb moisture from exhaled air and maintain its performance over longer periods of time. The advanced designs easily maintain their performance even during lengthy operations, while the less advanced models are intended for short periods of use only. Since used masks are considered heavily contaminated, it is imperative that</p> <ul style="list-style-type: none"> <li>- the wearer does not touch the main body of the mask with their fingers/hands,</li> <li>- hands are disinfected after removing the mask (complete hand disinfection),</li> <li>- the mask is worn over the nose and mouth; never wear the mask around the neck and</li> <li>- a used mask is disposed of if no longer required or between procedures; a new mask should be worn if further protection is required.</li> </ul>	Test	Result	Bacterial filtration efficiency (BFE), (%)	> 99%	Pressure difference (Pa/cm <sup>2</sup> )	< 60	Microbiological purity (KBE/g)	≤ 30	
Test	Result									
Bacterial filtration efficiency (BFE), (%)	> 99%									
Pressure difference (Pa/cm <sup>2</sup> )	< 60									
Microbiological purity (KBE/g)	≤ 30									

# Product data sheet

## wiroMEDI LOOP

VO (EU) 2023/988	Regulation (EU) No 2023/988 on general product safety
VO (EG) 1935/2004	Regulation regarding plastic materials and articles intended to come into contact with food.
VO (EU) 10/2011	Regulation on plastic materials and articles intended to come into contact with food

### Use

Instructions for use	<p><b>Before use</b></p> <p>Please ensure that you have selected the right product for your needs and areas of use. Do not wear disposable surgical face masks if they are faulty or have already been used. Disposable surgical face masks are intended for short-term use. In the medical sector, use is limited to one patient. Please follow the instructions for use in Annex A of the standard EN 14683:2019+AC:2019.</p> <p><b>Use</b></p> <p>Always wash your hands before putting on the mask. Do not touch your face while putting on the mask. Do not touch any part of the mask except the straps/ear loops and the outside of the mask. Avoid touching the inside of the mask. Masks should cover the nose and mouth completely and sit closely enough that as little air as possible can escape from the sides, top or bottom. Adjust the nosepiece to the shape of your nose.</p> <p>Please note that the mask is a standard size. Depending on the size and shape of your face or head, you may need to adjust the length of ear loops e.g. by knotting them, to make sure the mask fits properly.</p> <p>If you wear glasses, you can put them on again now. If your glasses fog up while wearing the mask, adjust the position of the mask on your face or treat your glasses with an anti-fogging product. Do not touch your mask while wearing it. Immediately replace a damaged or contaminated mask with an intact, clean one. Used masks are heavily contaminated and must be disposed of correctly. You should completely disinfect your hands after touching/removing the mask. Remove the mask if it is contaminated or damaged and dispose of it correctly.</p>
Safety notice	According to our risk analysis under EN ISO 14971:2019 + A11:2021, no additional safety instructions are required for this product.
Disposal	Unused and uncontaminated products can be incinerated in an environmentally friendly manner or disposed of at a waste disposal site. Contaminated products must be disposed of in accordance with the national laws and regulations valid for the respective type of contamination.
Shelf life	5 years from the production date, if unused and stored correctly.
Storage	Keep out of direct sunlight. Store in a dry place. Incorrect storage can lead to a reduction in the shelf life.
Traceability	Batch-related traceability from the determination of requirements to issuing. Product data sheet Version 2 Valid from Batch 2617 +

### Packaging

Recycling	Licensed in accordance with VerpackG (German Packaging Act), Licence number DE1618662887322. Please recycle the packaging.
-----------	---

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The numerical data for the measurements, weights, thicknesses and tolerances are target values, and there may be some deviation.

The information represents the current state of our knowledge. It offers no guarantee of the product characteristics and does not form the basis of a contractual legal understanding. Please reconfirm any information that is older than 2 years.

Typographical and other errors excepted.

The product data sheets provided are the property of WIROS Wilfried Rosbach GmbH. WIROS Wilfried Rosbach GmbH is the owner of all copyrights and rights of use and exploitation for the product data sheets. It is not permitted to change the product data sheets or to use the pictograms. The product data sheets provided may only be used in conjunction with original WIROS products.

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