[de]	Produktdatenblatt PP ASTROHAUBEN MIT MED. MASKE TYP II	Seite 2
[en]	Product data sheet PP ASTRO HOODS WITH MEDICAL FACE MASK TYPE II	Page 8
[nl]	Productgegevensformulier PP ASTROMUTSEN MET MED. MASKER TYPE II	Zijde 13
[fr]	Fiche produit CAGOULES ASTRO EN PP AVEC MASQUE MÉDICAL DE TYPE II	Page 19
[es]	Hoja de datos del producto GORROS INTEGRALES DE PP CON MASCARILLA MÉDICA TIPO II	Página 24
[it]	Scheda prodotto CAPPUCCI PROTETTIVI IN PP CON MASCHERINA MEDICA PROTETTIVA TIPO II	Pagina 30
[el]	Δελτίο δεδομένων προϊόντος ΣΚΟΥΦΟΣ PP ΤΥΠΟΥ ΑΣΤΡΟΝΑΥΤΗ ΜΕ ΙΑΤΡ. ΜΑΣΚΑ ΤΥΠΟΥ II	Σελίδα 36
[tr]	Ürün Bilgi Formu TIP II TIBBI MASKE İLE BİRLİKTE PP ASTRO BAŞLIK	Sayfa 42
[ru]	Техническая спецификация ШАПОЧКИ ТИПА ШЛЕМ С МАСКОЙ, МАТЕРИАЛ ПП МАСКА, ТИП II	Сторона 47
[pl]	Karta produktu KOMINIARKA TYPU ASTRO Z PP Z MEDYCZNĄ MASKĄ TYPU II	Strona 53
[da]	Produktdatablad PP ASTROHÆTTER MED MED. MUNDBIND TYPE II	Side 58
[sv]	Produktdatablad PP ASTROHUVA MED MEDICINSKT MUNSKYDD TYP II	Sida 63
[no]	Produktdatablad PP ASTROHETTER MED MED. MASKE TYPE II	Side 68
[cs]	Informační list výrobku PP ČEPICE ASTRO SE ZDRAV. MASKOU, TYP II	Strana 73
[hu]	Termékadatlap PP KÁMZSA II. TÍPUSÚ MASZK	Oldal: 78
[sk]	Informačný list výrobku PP KUKLY SO ZDRAV. MASKOU TYPU II	Strana 84
[sl]	Podatkovni list izdelka ASTRO KAPE PP Z MED. MASKO TIPA II	Stran 89
[ro]	Fișa produsului BONETĂ PP DE TIP ASTRONAUT, CU MASCĂ MED. TIP II	Pagina 94
[pt]	Folha de dados do produto TOUCA HOSPITALAR PP COM MÁSCARA MÉDICA TIPO II	Página 99
[lt]	Gaminio duomenų lapas PP APSAUGINIAI GALVOS GAUBTAI SU MED. KAUKĖ II TIPO	Puslapis 105
[lv]	Produkta datu lapa PP ASTRO KAPUCES AR MED. MASKU, II TIPS	Lappuse 110
[bg]	Информационен лист за продукта PP ASTRO ШАПКИ С МЕД. МАСКА ТИП II	Страница 115
[uk]	Технічний лист ПП ШАПОЧКИ-ШОЛОМ З МЕД. МАСКОЮ ТИПУ II	Сторінка 121
[et]	Toote teabeleht PP TUUKRIMÜTS KOOS MED. MASKIGA, TÜÜP II	Lk 126
[fi]	Tuotetietolehti PP ASTRONAUTIN HUPPU JA LÄÄK. MASKI TYYPPI II	Sivu 131

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PP ASTRO HOODS WITH MEDICAL FACE MASK TYPE II

Elastic neck cord | PP Medical face mask type II (3-ply | Elastic cord | Nosepiece)

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Authorisation

CH REP	Swiss AR Services AG, Industriestrasse 47, 6300 Zug, Switzerland
--------	--

Product description

Category	PREMIUM+
Product description	care & serve® PREMIUM+ wiroNOVA PP Astro hoods with medical face mask type II Elastic neck cord PP Medical face mask type II (3-ply Elastic cord Nosepiece)
Product features	3-ply nosepiece Elastic neck cord

Material

Material	Non-woven polypropylene
Material composition	Face mask: Outer layer: PP Spunbond
Material weight	Cap: ≈ 16 g/m ² (± 2 g/m ²) Face mask: Outer layer: ≈ 20 g/m ² (± 2 g/m ²) Filter medium: ≈ 22 g/m ² (± 2 g/m ²) Inner layer: ≈ 20 g/m ² (± 2 g/m ²)

Application

Intended use	Disposable cap to completely cover hair, nose and mouth to create a barrier.
Contact with food produce	Short-term indirect contact with foodstuffs is permitted.
Scope of application	There is no direct medical application on the patient. Medical protective clothing is intended to cover clothing and/or hair.
Application group	Employees and visitors from laboratories, the pharmaceutical industry, hospitals or similar institutions.
Patient group	Medical protective clothing is suitable for all patient groups.
Indication	Medical clothing as a barrier to cover clothing and/or hair.
Contraindication	No known contraindications.
Notes on use	Disposable product Not sterile
Limitation	The nose clip of the masks (face masks) is magnetic and not suitable in combination with highly magnetic devices, e.g. MRT devices. Not suitable for packing or storing food.
Note	All serious incidents directly related to this product must be reported to the manufacturer and the competent authorities.

Product data sheet

wiroNOVA

Application areas	Hospital Old people's home inc. specialist circles Medical practices Industry Pharmaceutical industry Laboratory Food processing	    
Supplementary material info	Latex-free	
Health advice	If skin irritation or allergic reactions occur, discontinue use immediately and seek medical advice. As far as is known according to the current state of the art, the products do not contain any toxic, carcinogenic or mutagenic substances when used in accordance with their intended use. These products do not contain latex components and are considered hypoallergenic. However, they can cause irritation in sensitive skin. Long-term wear may restrict skin ventilation, which can cause skin irritation.	

Standards

CE marking	Class I medical device in accordance with EU 2017/745	 								
SRN	DE-MF-000004956									
Basis UDI DI	4051642-wiroNOVA-D3									
	CE marking The CE mark guarantees the free movement of products and goods within the European Economic Community. All products with the CE mark comply with the essential requirements of the European regulations specified.									
EN 14683:2019+AC:2019	Medical face masks – requirements and test methods									
	<table border="1"> <thead> <tr> <th>Test</th> <th>Result</th> </tr> </thead> <tbody> <tr> <td>Bacterial filtration efficiency (BFE), (%)</td> <td>> 99 %</td> </tr> <tr> <td>Pressure difference (Pa/cm²)</td> <td>< 40</td> </tr> <tr> <td>Microbiological purity (KBE/g)</td> <td>≤ 30</td> </tr> </tbody> </table>	Test	Result	Bacterial filtration efficiency (BFE), (%)	> 99 %	Pressure difference (Pa/cm ²)	< 40	Microbiological purity (KBE/g)	≤ 30	
Test	Result									
Bacterial filtration efficiency (BFE), (%)	> 99 %									
Pressure difference (Pa/cm ²)	< 40									
Microbiological purity (KBE/g)	≤ 30									
	<p>Information</p> <p>User information</p> <p>When breathing, coughing and sneezing, smaller or larger quantities of respiratory droplet are released from the mucous membrane in the mouth and nose. Most droplets range from 0.5 µm to 12 µm in diameter. Larger droplets in particular can contain microorganisms from aerosolised particles. The droplets can then spread by air to sensitive places, such as an open surgical wound or sterile instruments.</p> <p>Medical masks, which are intended for use in operating theatres and healthcare facilities with similar requirements, are designed to protect the entire working environment. This standard describes two types of medical masks and their respective levels of protection. A Type I medical mask is used around patients to reduce the risk of spreading infection, particularly in epidemic or pandemic situations. A Type II mask is primarily intended for use by medical staff in an operating theatre or other medical facility with similar requirements.</p> <p>A special situation, which is also handled in the European Medical Device Regulation, is one in which the wearer requires spray protection against potentially contaminated liquids. If the mask is intended to protect the wearer against infectious germs (bacteria, viruses and fungi) when used correctly, respiratory protective equipment must be attached in accordance with Directive (89/686/EEC) and Regulation (EU) 2016/425 on personal protective equipment (PPE). The performance requirements of respiratory protective equipment are covered in EN 149.</p> <p>The effectiveness of a mask is dependent on a number of factors, such as filtration efficiency and how the mask sits on the wearer's face. Different designs suit different areas of application and so it is important to carefully choose a mask that produces the desired result.</p> <p>The filtration efficiency of mask materials can differ depending on the filter media. How the mask fits varies considerably between those who secure it behind their ears using ear loops and those who secure it using straps over the head and a nose piece that can be adjusted to fit the shape of the face. Mask performance is typically determined by in-vitro tests of the materials from which the mask is made. However, it is important to carefully consider how the mask fits when choosing a mask for a specific purpose. Another factor to consider is the mask's ability to absorb moisture from exhaled air and maintain its performance over longer periods of time. The advanced designs easily maintain their performance even during lengthy operations, while the less advanced models are intended for short periods of use only. Since used masks are considered heavily contaminated, it is imperative that</p> <ul style="list-style-type: none"> - the wearer does not touch the main body of the mask with their fingers/hands, - hands are disinfected after removing the mask (complete hand disinfection), - the mask is worn over the nose and mouth; never wear the mask around the neck and - a used mask is disposed of if no longer required or between procedures; a new mask should be worn if further protection is required. 									
VO (EU) 2023/988	Regulation (EU) No 2023/988 on general product safety									

Product data sheet

wiroNOVA

VO (EG) 1935/2004 Regulation regarding plastic materials and articles intended to come into contact with food.

VO (EU) 10/2011 Regulation on plastic materials and articles intended to come into contact with food

Use

Instructions for use

Before use
Make sure that you have selected the right product for your needs and areas of use. Do not wear items of protective medical clothing if they are faulty or have already been used. Items of disposable protective medical clothing are intended for short-term use. In the medical sector, use is limited to one patient.

Handling the hood
Pull the hood over your head. Make sure that your hair is completely covered by the hood. This also applies to long hair.

Handling the face mask
Please follow the instructions for use in attachment A EN 14683:2019+AC:2019.
Always wash your hands before putting on the mask. Do not touch your face while putting on the mask. Do not touch any part of the mask except the straps/ear loops and the outside of the mask. Do not touch the inside of the mask. Masks should cover the nose and mouth and sit closely enough that as little air as possible escapes from the sides, top or bottom. Adjust the nosepiece around your nose.
Please note that the mask is a standard size. Depending on the size and shape of your face or head, you may need to adjust the length of ear loops e.g. by knotting them, to make sure the mask fits properly.
If you wear glasses, you can now put them on again. If your glasses fog up while wearing the mask, adjust the position of the mask on your face or treat your glasses with an anti-fogging product. Do not touch your mask while wearing it. Immediately replace a damaged or contaminated mask with an intact, clean mask. Used masks are heavily contaminated and must be disposed of correctly. You should completely disinfect your hands after touching/removing the mask. Remove the mask if it is contaminated or damaged and dispose of it correctly.

Safety notice According to our risk analysis under EN ISO 14971:2019 + A11:2021, no additional safety instructions are required for this product.

Disposal Unused and uncontaminated products can be incinerated in an environmentally friendly manner or disposed of at a waste disposal site. Contaminated products must be disposed of in accordance with the national laws and regulations valid for the respective type of contamination.

Shelf life 5 years from the production date, if unused and stored correctly.

Storage Keep out of direct sunlight.
Store in a dry place.
Incorrect storage can lead to a reduction in the shelf life.

Traceability Batch-related traceability from the determination of requirements to issuing.
Product data sheet Version 2 Valid from Batch 2617 +

Packaging

Recycling Licensed in accordance with VerpackG (German Packaging Act), Licence number DE1618662887322.
Please recycle the packaging.

www.wiros.de

WEB <https://www.wiros.de/shop/kopfbedeckungen/pp-astrohauben/wironova/>

Variant 1

Sales unit 1000 Unit/Cardboard box

Packaging level 1 100 Unit/Polybag

Minimum order quantity 1 Cardboard box

Product data sheet

wiroNOVA

Item numbers	Item number	Colour	Size			
	10002509	white	Unisize			
	10002519	green	Unisize			
	10002529	blue	Unisize			
	10002549	yellow	Unisize			

GTINs	Item number	Unit of Use UDI	Packaging level 1 UDI	Packaging level 2 UDI	Cardboard box UDI	
	10002509	4051642009162	4051642009155	-----	4051642009148	
	10002519	4051642009193	4051642009186	-----	4051642009179	
	10002529	4051642009223	4051642009216	-----	4051642009209	
	10002549	4051642009254	4051642009247	-----	4051642009230	

The numerical data for the measurements, weights, thicknesses and tolerances are target values, and there may be some deviation.

The information represents the current state of our knowledge. It offers no guarantee of the product characteristics and does not form the basis of a contractual legal understanding. Please reconfirm any information that is older than 2 years.

Typographical and other errors excepted.

The product data sheets provided are the property of WIROS Wilfried Rosbach GmbH. WIROS Wilfried Rosbach GmbH is the owner of all copyrights and rights of use and exploitation for the product data sheets. It is not permitted to change the product data sheets or to use the pictograms. The product data sheets provided may only be used in conjunction with original WIROS products.

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