[de]	Produktdatenblatt ZAHNBÜRSTEN	Seite 2
[en]	Product data sheet TOOTHBRUSHES	Page 5
[nl]	Productgegevensformulier TANDENBORSTELS	Zijde 8
[fr]	Fiche produit BROSSES À DENTS	Page 11
[es]	Hoja de datos del producto CEPILLOS DE DIENTES	Página 14
[it]	Scheda prodotto SPAZZOLINI DA DENTI	Pagina 17
[el]	Δελτίο δεδομένων προϊόντος ΟΔΟΝΤΟΒΟΥΡΤΣΕΣ	Σελίδα 20
[tr]	Ürün Bilgi Formu DIŞ FIRÇALARI	Sayfa 23
[ru]	Техническая спецификация ЗУБНЫЕ ЩЕТКИ	Сторона 26
[pl]	Karta produktu SZCZOTECZKI DO ZĘBÓW	Strona 29
[da]	Produktdatablad TANDBØRSTER	Side 32
[sv]	Produktdatablad TANDBORSTAR	Sida 35
[no]	Produktdatablad TANNBØRSTER	Side 38
[cs]	Informační list výrobku ZUBNÍ KARTÁČKY	Strana 41
[hu]	Termékadatlap FOGKEFÉK	Oldal: 44
[sk]	Informačný list výrobku ZUBNÉ KEFKY	Strana 47
[sl]	Podatkovni list izdelka ZOBNE ŠČETKE	Stran 50
[ro]	Fișa produsului PERIUȚE DE DINȚI	Pagina 53
[pt]	Folha de dados do produto ESCOVAS DE DENTES	Página 56
[lt]	Gaminio duomenų lapas DANTŲ ŠEPETĖLIAI	Puslapis 59
[lv]	Produkta datu lapa ZOBU SUKAS	Lappuse 62
[bg]	Информационен лист за продукта ЧЕТКИ ЗА ЗЪБИ	Страница 65
[uk]	Технічний лист ЗУБНІ ЩІТКИ	Сторінка 68
[et]	Toote teabeleht HAMBAHARJAD	Lk 71
[fi]	Tuotetietolehti HAMMASHARJAT	Sivu 74


ZAHNBÜRSTEN

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PREMIUM+

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wiroDENTAL

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
Product data sheet

wiroDENTAL

Manufacturer	
Company + Address	WIROS Wilfried Rosbach GmbH, Karl-Arnold-Str. 5, 47877 Willich, Germany
Trademark	care & serve®
Product description	
Category	PREMIUM+
Product description	care & serve® PREMIUM+ wiroDENTAL Toothbrushes Integrated toothpaste
Product features	Integrated toothpaste
Application	
Intended use	Disposable toothbrushes for brushing teeth without the need for water or toothpaste.
Notes on use	Disposable product 
Application areas	Hospital Old people's home inc. specialist circles Medical practices Cosmetics Wellness Beauty industry   
Supplementary material info	Latex-free 
Standards	
VO (EU) 2023/988	Regulation (EU) No 2023/988 on general product safety
Use	
Instructions for use	Before use Remove the disposable toothbrush from the packaging and check it for damage or defects. Make sure that it is clean and hygienic. Use Moisten the bristles on the disposable toothbrush under running water, if possible, to make brushing easier. As the toothpaste is integrated into the bristles, no additional toothpaste is needed. Move the disposable toothbrush at a 45 degree angle to your teeth. Brush your teeth in gentle, circular motions or with short up and down movements. Brush all the tooth surfaces, including the outer, inner and chewing surfaces as well as your tongue and gums. Brush your teeth for at least two minutes to make sure that all areas are cleaned thoroughly. Rinse your mouth out thoroughly with water in order to remove any food residue and toothpaste. Disposable toothbrushes are designed for single use and should be disposed of after use.
Safety notice	
Disposal	Unused and uncontaminated products can be incinerated in an environmentally friendly manner or disposed of at a waste disposal site. Contaminated products must be disposed of in accordance with the national laws and regulations valid for the respective type of contamination.
Shelf life	5 years from the production date, if unused and stored correctly.
Storage	Keep out of direct sunlight. Store in a dry place. Incorrect storage can lead to a reduction in the shelf life.  
Traceability	Batch-related traceability from the determination of requirements to issuing. Product data sheet Version 2 Valid from Batch 2617 +

Product data sheet

wiroDENTAL

Packaging

Recycling	Licensed in accordance with VerpackG (German Packaging Act), Licence number DE1618662887322. Please recycle the packaging.	 VerpackG
-----------	---	--

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
The numerical data for the measurements, weights, thicknesses and tolerances are target values, and there may be some deviation.
The information represents the current state of our knowledge. It offers no guarantee of the product characteristics and does not form the basis of a contractual legal understanding. Please reconfirm any information that is older than 2 years.
Typographical and other errors excepted.
The product data sheets provided are the property of WIROS Wilfried Rosbach GmbH. WIROS Wilfried Rosbach GmbH is the owner of all copyrights and rights of use and exploitation for the product data sheets. It is not permitted to change the product data sheets or to use the pictograms. The product data sheets provided may only be used in conjunction with original WIROS products.

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